

## VELP Scientifica Solutions for Stirring



**MST**  
(Analog / Digital)

**MICROSTIRRER**

**ESP**

**MULTI-HS 6**

**MULTI-HS 15**

**MULTISTIRRER 6**  
(Analog / Digital)

**MULTISTIRRER 15**  
(Analog / Digital)

**AMI**

**AMI4**

Stirring speed rpm	up to 1100	up to 1100	up to 1100	up to 1000	up to 1000	from 80 to 1500	from 80 to 1500	up to 1100	up to 1100	
Stirring volume L	up to 5	up to 5	up to 5	6 x 100ml	15 x 74 ml	up to 2.4 (400 ml each)	up to 3.75 (250 ml each)	up to 5	up to 20	
Timer	• (digital version only)			Yes	Yes	• (digital version only)	• (digital version only)			
Reverse Operation				Yes	Yes	• (digital version only)	• (digital version only)			
Dimensions (WxHxD) mm	120x50x145	120x48x128	160x33x230	240x65x550	240x65x550	230x51.5x370	230x51.5x370	150x55x270	600x55x270	
(in)	(4.7x1.97x5.9)	(4.7x1.9x5.0)	(6.3x1.3x9.0)	(9.4x2.6x21.7)	(9.4x2.6x21.7)	(9.0x2.0x14.5)	(9.0x2.0x14.5)	(5.9x2.2x10.6)	(23.6x2.2x10.6)	
Weight	Kg	0.6	0.55	0.9	4,7	4,8	1.75	2.1	1.2	4.0
(lb)	(1.3)	(1.1)	(2.0)	(10,4)	(10,5)	(3.8)	(4.6)	(2.6)	(8.8)	
Power Supply	100÷240 V	100÷240 V	100÷240 V	115 or 230 V	115 or 230 V	100÷240 V	100÷240 V	100÷240 V	100÷240 V	
Power	0.6 W	0.6 W	5 W	600 W	600 W	3.6 W	9 W	1.2 W	4.8 W	
Temperature	-	-	-	Up to 120°C	Up to 120°C	-	-	-	-	



www.jointlab.com  
www.frigolab.it  
Via C. Treves 57, Trezzano sul Naviglio  
Milano, Italy  
+39.02.39310823  
info@jointlab.com