

## VELP Scientifica MSL Digital Series



MSL 8 Digital



MSL 25 Digital



MSL 50 Digital

|                         |            |            |            |
|-------------------------|------------|------------|------------|
| Stirring speed rpm      | 40 to 1500 | 40 to 1500 | 20 to 1100 |
| Stirring volume L (H2O) | up to 8    | up to 25   | up to 50   |
| Timer                   | •          | •          | •          |
| Reverse Operation       | •          | •          | •          |
| Intermittent Mode       | •          | •          | •          |
| Dimensions (WxHxD) mm   | 168x60x183 | 253x71x269 | 354x77x369 |
| Weight Kg               | 0.95       | 1.8        | 4          |
| Power Supply            | 100÷240 V  | 100÷240 V  | 100÷240 V  |
| Power input             | 5 W        | 9 W        | 10 W       |



[www.jointlab.com](http://www.jointlab.com)

[www.frigolab.it](http://www.frigolab.it)

Via C. Treves 57, Trezzano sul Naviglio

Milano, Italy

+39.02.39310823

[info@jointlab.com](mailto:info@jointlab.com)