


RX3



ZX3



ZX4



TX4



CLASSIC



WIZARD

Stirring speed rpm		3000, constant	0 ÷ 3000	0 ÷ 3000	0 ÷ 3000	0 ÷ 3000	0 ÷ 3000
Orbital diameter mm		4.5	4.5	4.5	4.5	4.5	4.5
Operating mode	touch	●	●			●	
	continuous		●	●	●	●	●
	IR sensor			●	●		●
Digital display and timer					●		
Dimensions (WxHxD) mm		150x130x165	150x130x165	150x130x165	150x130x165	180x70x220	180x70x220
	(in)	(5.9x5.1x6.5)	(5.9x5.1x6.5)	(5.9x5.1x6.5)	(5.9x5.1x6.5)	(7.1x2.8x8.7)	(7.1x2.8x8.7)
Weight	Kg	2.7	2.7	2.7	2.7	2.2	2.2
	(lb)	(5.9)	(5.9)	(5.9)	(5.9)	(4.9)	(4.9)
Power Supply		100 ÷ 240 V	100 ÷ 240 V	100 ÷ 240 V	100 ÷ 240 V	100 ÷ 240 V	100 ÷ 240 V
Power		15 W	15 W	15 W	15 W	15 W	15 W